


NASAL-BILE-DUCT DRAINAGE-SET, TYPE NAG

with 0.035 " guide wire, 7fr. nasal
tube connector, funnel 3-way stopcock



SKU: PK-KNB7NAG
Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Single Use Instruments](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Single Use Instruments](#)