

NASAL-BILE-DUCT DRAINAGE-SET, TYPE PARIS

with 0.035 " guide wire, 10 fr. nasal

tube connector, funnel 3-way stopcock



SKU: PK-KNBL10
Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Single Use Instruments](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Single Use Instruments](#)