

## NEEDLE CASE

perforated

65x42x8.0mm



**SKU:** 250-304-65

**Categories:** [Open Surgery](#), [Various Items](#)

## ZUSÄTZLICHE INFORMATIONEN

|                   |                       |
|-------------------|-----------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>  |
| <b>TYPE</b>       | <a href="#">Bowls</a> |
| <b>KLS Martin</b> | 56-302-65             |

## ZUSÄTZLICHE INFORMATIONEN

|                   |                       |
|-------------------|-----------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>  |
| <b>TYPE</b>       | <a href="#">Bowls</a> |
| <b>KLS Martin</b> | 56-302-65             |