

NEEDLE CASE

perforated

75x25x4.5mm



SKU: 250-304-75

Categories: [Open Surgery](#), [Various Items](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------|
| UNIT | pcs. |
| TYPE | Bowls |
| Aesculap | BL915 |
| KLS Martin | 56-302-75 |
| Medicon | 46.37.48 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------|
| UNIT | pcs. |
| TYPE | Bowls |
| Aesculap | BL915 |
| KLS Martin | 56-302-75 |
| Medicon | 46.37.48 |