

## NEEDLE CASE

perforated

75x25x4.5mm



**SKU:** 250-304-75

**Categories:** [Open Surgery](#), [Various Items](#)

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bowls</a>
<b>Aesculap</b>	BL915
<b>KLS Martin</b>	56-302-75
<b>Medicon</b>	46.37.48

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bowls</a>
<b>Aesculap</b>	BL915
<b>KLS Martin</b>	56-302-75
<b>Medicon</b>	46.37.48