

NEEDLE ELECTRODE; FLEXIBLE

bipolar; 5 Fr.; 360mm

needle angled

SKU: 582-516-05

Categories: [Endoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#), [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [HF-Electrodes](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [HF-Electrodes](#)