

NEEDLE ELECTRODE; FLEXIBLE

monopolar; 5 Fr.; 470mm

SKU: 582-511-05

Categories: [Endoscopy](#), [GUBBINI/Gyn](#), [Gynecology](#), [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [HF-Electrodes](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [HF-Electrodes](#)