

## NEEDLE HOLDER Ø 5MM, CONVEX/ CONCAVE

with Aluminium handle, teflon coated

WL 450mm, LL-Connection



**SKU:** 614-120-54

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#)

**TYPE** [LOFTLINE](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Needle Holders](#)

**TYPE** [LOFTLINE](#)