

NEEDLE RECEIVING FORCEPS; Ø 5MM

with Aluminium handle, teflon coated

WL 330mm, LL-Connection



SKU: 614-195-53
Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Needle Holders
TYPE	LOFTLINE

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Needle Holders
TYPE	LOFTLINE