

## NEIVERT NEEDLE HOLDER

straight; smooth; special groove

13.0cm (5.2")




**SKU:** 050-046-13

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>KLS Martin</b>	20-338-13
<b>AUTHORS</b>	<a href="#">NEIVERT</a>
<b>Medicon</b>	10.08.12
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>KLS Martin</b>	20-338-13
<b>AUTHORS</b>	<a href="#">NEIVERT</a>
<b>Medicon</b>	10.08.12

## SPECIFICATIONS

[straight](#)