

NEPHROSCOPE; 45° ANGELD

view 15°;working channel Ø 5.0mm

15 Charr.; working lenght 250mm

ACMI- Wolf- and Storz adapters

SKU: 601-140-45

Categories: [Endoscopy](#), [Equipment](#), [Urology](#)




ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Telescopes](#)

SCOPE [Urology](#)