

NEPHROSCOPE SHEATH 26 FR.

FL 264mm; 1 x turnable stopcock

colour code: black



SKU: 630-420-44

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|-----------------------------|
| UNIT | pcs. |
| RANGE | Nephroscopy |
| TYPE | Sheaths |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|-----------------------------|
| UNIT | pcs. |
| RANGE | Nephroscopy |
| TYPE | Sheaths |