

## NITINOL-BASKET , TIPLESS ROUND

hexagonal; 30 x 60mm; 4 wires

2.6mm; 230cm



**SKU:** PK-1NRPK2630230K-Z

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)