

OBTURATOR VISUAL FOR RESECTOSCOPY

24 / 26Fr.; with 5Fr.working channel
for 4mm/30° telescope



SKU: 630-230-24
Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Obturators

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Obturators