

OBTURATOR VISUAL FOR RESECTOSCOPY

27 / 28.5Fr.; with 7Fr. working channel



SKU: 630-230-27
Categories: [Endoscopy](#), [Gynecology](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Obturators

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Resectoscopy 4.0
TYPE	Obturators