

OLIVECRONA TRIGEMINAL SCISSORS

serrated

working length: 120 mm




SKU: 010-469-23

Categories: [Open Surgery](#), [Scissors](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Scissors for Neuro Surgery , Scissors for Vascular Surgery |
| Aesculap | BC665 |
| KLS Martin | 11-475-23 |
| AUTHORS | OLIVECRONA |
| Medicon | 57.25.70 |
| SPECIFICATIONS | straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Scissors for Neuro Surgery , Scissors for Vascular Surgery |
| Aesculap | BC665 |
| KLS Martin | 11-475-23 |

AUTHORS

[OLIVECRONA](#)

Medicon

57.25.70

SPECIFICATIONS

[straight](#)