

OLSEN-HEGAR NEEDLE HOLDER

straight; serrated; with scissors

17.0cm (6.8")




SKU: 050-310-17

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM291
KLS Martin	20-460-17
AUTHORS	HEGAR , OLSEN
Medicon	10.60.17
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM291
KLS Martin	20-460-17

AUTHORS

[HEGAR](#), [OLSEN](#)

Medicon

10.60.17

SPECIFICATIONS

[serrated](#), [straight](#)