

## OTOLOGY SUCTION TUBE

tapered; curved; Luer connector

Ø 3.5mm; WL: 11.5cm

**SKU:** 086-500-35

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Suction Tubes](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Suction Tubes](#)