


OVERHOLT LIGATURE FORCEPS

curved; serrated

29.5cm (11.5")



SKU: 040-427-01
 Categories: [Open Surgery](#), [Hemostats](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ091
KLS Martin	13-449-01
AUTHORS	OVERHOLT
SPECIFICATIONS	curved , serrated

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dissecting and Ligature Forceps
Aesculap	BJ091
KLS Martin	13-449-01
AUTHORS	OVERHOLT

SPECIFICATIONS

[curved](#), [serrated](#)