

## PENNINGTON HAEMOSTATIC FORCEPS

straight; serrated

15.5cm (6.2")




**SKU:** 040-358-15

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH902
<b>KLS Martin</b>	13-280-15
<b>AUTHORS</b>	<a href="#">PENNINGTON</a>
<b>Medicon</b>	17.34.15
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a>
<b>Aesculap</b>	BH902
<b>KLS Martin</b>	13-280-15

## AUTHORS

[PENNINGTON](#)

## Medicon

17.34.15

## SPECIFICATIONS

[serrated](#), [straight](#)