

PHANEUF HYSTERECTOMY FORCEPS

straight; 1X2 teeth

21.5cm (8.6")



SKU: 040-630-21

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Hysterectomy Forceps |
| Aesculap | BJ500 |
| KLS Martin | 13-600-21 |
| AUTHORS | PHANEUF |
| Medicon | 16.80.21 |
| SPECIFICATIONS | straight , toothed |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------------|
| UNIT | pcs. |
| TYPE | Hysterectomy Forceps |
| Aesculap | BJ500 |
| KLS Martin | 13-600-21 |

AUTHORS

[PHANEUF](#)

Medicon

16.80.21

SPECIFICATIONS

[straight](#), [toothed](#)