

## PLESTER KNIFE

Ø 3.5mm

16.0cm (6.4")

**SKU:** 111-671-16

**Categories:** [Open Surgery](#), [Otology](#)



## ZUSÄTZLICHE INFORMATIONEN

|                   |                         |
|-------------------|-------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>    |
| <b>TYPE</b>       | <a href="#">Knives</a>  |
| <b>Aesculap</b>   | OG106                   |
| <b>KLS Martin</b> | 36-020-21               |
| <b>AUTHORS</b>    | <a href="#">PLESTER</a> |



## ZUSÄTZLICHE INFORMATIONEN

|                   |                         |
|-------------------|-------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>    |
| <b>TYPE</b>       | <a href="#">Knives</a>  |
| <b>Aesculap</b>   | OG106                   |
| <b>KLS Martin</b> | 36-020-21               |
| <b>AUTHORS</b>    | <a href="#">PLESTER</a> |