

PLESTER SELFRETAINING RETRACTOR

sharp; 2X2 teeth
 13.0cm (5.2")



SKU: 070-830-13
Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV078
KLS Martin	15-756-13
AUTHORS	PLESTER
SPECIFICATIONS	sharp , toothed

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV078
KLS Martin	15-756-13
AUTHORS	PLESTER

SPECIFICATIONS

[sharp](#), [toothed](#)