

PLESTER SELFRETAINING RETRACTOR

blunt; 2X2 teeth
 13.0cm (5.2")



SKU: 070-831-13
Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV079
KLS Martin	15-755-13
AUTHORS	PLESTER
Medicon	24.26.03
SPECIFICATIONS	blunt , toothed

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV079
KLS Martin	15-755-13

AUTHORS

[PLESTER](#)

Medicon

24.26.03

SPECIFICATIONS

[blunt](#), [toothed](#)