

PLESTER SUCTION TUBE

7 Fr.; 30° angled

14.0cm (5.2") working length

SKU: 085-032-07

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| AUTHORS | PLESTER |

ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| AUTHORS | PLESTER |