

## POLITZER TYMPANUM NEEDLE

angled; needle vertically;

16.5cm (6.6")



**SKU:** 111-616-16

**Categories:** [Open Surgery](#), [Otology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Paracentesis Needles, Knives</a>
<b>Aesculap</b>	OF608
<b>KLS Martin</b>	36-655-16
<b>AUTHORS</b>	<a href="#">POLITZER</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Paracentesis Needles, Knives</a>
<b>Aesculap</b>	OF608
<b>KLS Martin</b>	36-655-16
<b>AUTHORS</b>	<a href="#">POLITZER</a>