

## POLYPECTOMY LOOP

hexagonal 20mm

2.3mm, length 230cm



**SKU:** PK-4PHPK2320230

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Polypectomy Loops, EMR-Snares](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Polypectomy Loops, EMR-Snares](#)