


POLYPECTOMY LOOP

oval; 30mm, spring flat reinforced

tubus (7fr.) 2.3 x 250cm; single use



SKU: PK-4EOPK2330250F-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)