


## POLYPECTOMY LOOP, ISOLATED TIP

oval 10mm

2.3mm, length: 230cm



**SKU:** PK-4POPK2310230ID

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Polypectomy Loops, EMR-Snares](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Polypectomy Loops, EMR-Snares](#)