


## POLYPECTOMY LOOP , MONOFIL (SINGLE USE)

snare: oval; 15mm, box 5/pcs

tubus 2.3 x 230cm



**SKU:** PK-4EOPK2315230M-S

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Single Use Instruments](#)