


POLYPECTOMY LOOP

hexagonal; 25mm

tubus (5fr.) 1.8x250cm; single use



SKU: PK-4EHPK1825250-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)