

POOLE SUCTION TUBE; STRAIGHT

perforated and removable external tube

inner Ø 4.0mm; outer Ø 8.0mm; 24.5cm

SKU: 085-020-08

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	POOLE

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	POOLE