

## POOLE SUCTION TUBE; STRAIGHT

perforated and removable external tube

inner Ø 7.0mm; outer Ø 10.0mm; 23.5cm

**SKU:** 085-020-10

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>Aesculap</b>	GF862
<b>KLS Martin</b>	18-550-22
<b>AUTHORS</b>	<a href="#">POOLE</a>
<b>Medicon</b>	09.87.30

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>Aesculap</b>	GF862
<b>KLS Martin</b>	18-550-22
<b>AUTHORS</b>	<a href="#">POOLE</a>

**Medicon**

09.87.30