

POOLE SUCTION TUBE; STRAIGHT

perforated and removable external tube

inner Ø 3.7mm; outer Ø 6.0mm; 19.0cm

SKU: 085-020-06

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	POOLE
Medicon	09.87.18

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
AUTHORS	POOLE
Medicon	09.87.18