

POWERGRIP BIPOLAR COAGULATION FORCEPS

ring handle with flat plug
(ex.628-301-53)



SKU: 628-302-52
Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Handles

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Handles