

## PRESBYTERIAN TUBING CLAMP

straight; smooth

18.0cm (7.2")




**SKU:** 032-012-18

**Categories:** [Open Surgery](#), [Hemostats](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>AUTHORS</b>	<a href="#">PRESBYTERIAN</a>
<b>Medicon</b>	17.72.18
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>AUTHORS</b>	<a href="#">PRESBYTERIAN</a>
<b>Medicon</b>	17.72.18
<b>SPECIFICATIONS</b>	<a href="#">straight</a>