

## PRESBYTERIAN TUBING CLAMP

straight; smooth

20.0cm (8")




**SKU:** 032-012-20

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>AUTHORS</b>	<a href="#">PRESBYTERIAN</a>
<b>Medicon</b>	17.72.20
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>AUTHORS</b>	<a href="#">PRESBYTERIAN</a>
<b>Medicon</b>	17.72.20
<b>SPECIFICATIONS</b>	<a href="#">straight</a>