

## PRESBYTERIAN TUBING CLAMP

straight; smooth

16.5cm (6.6")




**SKU:** 032-012-16

**Categories:** [Open Surgery](#), [Hemostats](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>Aesculap</b>	AN909
<b>KLS Martin</b>	13-942-16
<b>AUTHORS</b>	<a href="#">PRESBYTERIAN</a>
<b>Medicon</b>	17.72.16
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>Aesculap</b>	AN909
<b>KLS Martin</b>	13-942-16

## AUTHORS

[PRESBYTERIAN](#)

## Medicon

17.72.16

## SPECIFICATIONS

[straight](#)