

PRESBYTERIAN TUBING CLAMP

straight; smooth

16.5cm (6.6")




SKU: 032-012-16

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
Aesculap	AN909
KLS Martin	13-942-16
AUTHORS	PRESBYTERIAN
Medicon	17.72.16
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
Aesculap	AN909
KLS Martin	13-942-16

AUTHORS

[PRESBYTERIAN](#)

Medicon

17.72.16

SPECIFICATIONS

[straight](#)