


## PROBE

Ø 2.0mm

11.5cm (4.6")



**SKU:** 080-004-11

**Categories:** [Open Surgery](#), [Probes](#)

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Buttoned Probes</a>
<b>KLS Martin</b>	16-100-11
<b>Medicon</b>	09.05.11
<b>SPECIFICATIONS</b>	<a href="#">without eye</a>

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Buttoned Probes</a>
<b>KLS Martin</b>	16-100-11
<b>Medicon</b>	09.05.11
<b>SPECIFICATIONS</b>	<a href="#">without eye</a>