


PROBE

Ø 2.0mm

14.5cm (5.8")



SKU: 080-004-14
Categories: [Open Surgery](#), [Probes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN134
KLS Martin	16-100-14
Medicon	09.05.14
SPECIFICATIONS	without eye

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN134
KLS Martin	16-100-14
Medicon	09.05.14

SPECIFICATIONS

[without eye](#)