


PROBE

Ø 2.0mm

18.0cm (7.2")



SKU: 080-004-18
Categories: [Open Surgery](#), [Probes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN138
KLS Martin	16-100-18
Medicon	09.05.18
SPECIFICATIONS	without eye

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN138
KLS Martin	16-100-18
Medicon	09.05.18

SPECIFICATIONS

[without eye](#)