


PROBE

Ø 2.0mm

20.0cm (8")



SKU: 080-004-20

Categories: [Open Surgery](#), [Probes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN140
KLS Martin	16-100-20
Medicon	09.05.20
SPECIFICATIONS	without eye

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN140
KLS Martin	16-100-20
Medicon	09.05.20

SPECIFICATIONS

[without eye](#)