


PROBE

Ø 1.0mm

14.5cm (5.8")



SKU: 080-000-14

Categories: [Open Surgery](#), [Probes](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN104
KLS Martin	16-102-14
Medicon	09.03.14
SPECIFICATIONS	without eye

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN104
KLS Martin	16-102-14
Medicon	09.03.14

SPECIFICATIONS

[without eye](#)