


## PROBE

Ø 1.0mm

14.5cm (5.8")



**SKU:** 080-000-14

**Categories:** [Open Surgery](#), [Probes](#)

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Buttoned Probes</a>
<b>Aesculap</b>	BN104
<b>KLS Martin</b>	16-102-14
<b>Medicon</b>	09.03.14
<b>SPECIFICATIONS</b>	<a href="#">without eye</a>

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Buttoned Probes</a>
<b>Aesculap</b>	BN104
<b>KLS Martin</b>	16-102-14
<b>Medicon</b>	09.03.14

## SPECIFICATIONS

[without eye](#)