

PROBE WITH EYE

Ø 2.0mm

14.5cm (5.8")

SKU: 080-013-14

Categories: [Open Surgery](#), [Probes](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN034
KLS Martin	16-108-14
Medicon	09.10.14
SPECIFICATIONS	with eye

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Buttoned Probes
Aesculap	BN034
KLS Martin	16-108-14
Medicon	09.10.14

SPECIFICATIONS

[with_eye](#)