

PROBE WITH EYE

Ø 2.0mm

18.0cm (7.2")

SKU: 080-013-18

Categories: [Open Surgery](#), [Probes](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Buttoned Probes](#)

KLS Martin 16-108-18

SPECIFICATIONS [with eye](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Buttoned Probes](#)

KLS Martin 16-108-18

SPECIFICATIONS [with eye](#)