

## PROBE WITH EYE

Ø 2.0mm

20.0cm (8")

**SKU:** 080-013-20

**Categories:** [Open Surgery](#), [Probes](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Buttomed Probes</a>
<b>Aesculap</b>	BN040
<b>KLS Martin</b>	16-108-20
<b>Medicon</b>	09.10.20
<b>SPECIFICATIONS</b>	<a href="#">with eye</a>

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Buttomed Probes</a>
<b>Aesculap</b>	BN040
<b>KLS Martin</b>	16-108-20
<b>Medicon</b>	09.10.20

## SPECIFICATIONS

[with\\_eye](#)