

## QUINKE LUMBAR PUNCTURE NEEDLE

Ø 0.7 x 80mm; Luer-Lock connection

Minimum order quantity: 10 pieces

**SKU:** 086-120-05

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Biopsy Needles</a>
<b>KLS Martin</b>	18-630-18
<b>AUTHORS</b>	<a href="#">QUINKE</a>
<b>Medicon</b>	09.72.07

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Biopsy Needles</a>
<b>KLS Martin</b>	18-630-18
<b>AUTHORS</b>	<a href="#">QUINKE</a>
<b>Medicon</b>	09.72.07