

## RESECTOSCOPE SHEATH 9 FR.

with rotatable stopcock

distal insulated, colour code blue



**SKU:** 630-512-13

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Pediatric Resectoscopy</a>
<b>TYPE</b>	<a href="#">Sheaths</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Pediatric Resectoscopy</a>
<b>TYPE</b>	<a href="#">Sheaths</a>